Bucknell University Bucknell Digital Commons

Other Faculty Research and Publications

Faculty Scholarship

Fall 9-2019

Technical Report 2019-01: Pupil Labs Eye Tracking User Guide

Joan D. Gannon Bucknell University

Augustine Ubah Bucknell University, acu002@bucknell.edu

Chris Dancy Bucknell University, cld028@bucknell.edu

Follow this and additional works at: https://digitalcommons.bucknell.edu/fac_pubs

Part of the Artificial Intelligence and Robotics Commons, Cognitive Psychology Commons, and the Other Computer Sciences Commons

Recommended Citation

Gannon, J., Ubah, A., & Dancy, C. L. (2019). Technical Report 2019-01: Pupil Labs Eye Tracking User Guide. Bucknell University Department of Computer Science Technical Report Series. Computer Science. Bucknell University. Lewisburg, PA.

This Unpublished Paper is brought to you for free and open access by the Faculty Scholarship at Bucknell Digital Commons. It has been accepted for inclusion in Other Faculty Research and Publications by an authorized administrator of Bucknell Digital Commons. For more information, please contact dcadmin@bucknell.edu.

Department of Computer Science Bucknell University

Technical Report 2019-01

Pupil Labs Eye Tracking User Guide

Joanie Gannon, Augustine Ubah, Christopher L. Dancy Bucknell University

September 07, 2019

Table of Contents

Getting set up	4
The Hardware	4
The Software	4
Using the software	6
Pupil Service	6
How to navigate this app	7
Pupil Capture	10
Pupil Player	16
Getting good data	. 19
Troubleshooting	. 22
I'm getting frozen images/a grey screen in Pupil Service	22
Pupil Capture's calibration is not working	22
Links and References	. 24

Introduction

The purpose of this document is to allow users to quickly set up and capture data from the pupil labs eye tracking glasses without any previous experience.



For additional resources, see the <u>Links and References</u> section of this guide.

Getting set up

To get started, you will need to set up the physical hardware for the glasses, as well as install the necessary software. More information on setup process can be found on this website: https://docs.pupil-labs.com/#getting-started

The Hardware

Included with the glasses, you should have the following: -A USB cable -2 Orange extenders -An additional nose pads -An extra camera lens

The Pupil Labs glasses are relatively easy to set up. Simply connect one end of the USB cable to the computer, and the other end to the glasses.



More information about how to adjust the glasses to properly capture the eye is in the <u>Getting Good Data</u> section. It is much easier to perform these adjustments once the software is installed.

The Software

To begin, download the companion software by clicking the following link:

https://github.com/pupil-labs/pupil/releases/latest (At time of writing, the latest version of this software is 1.8) A screen similar to the one shown in the image below should show up, and select the appropriate download for your operating system.



After the file is downloaded, unzip it and open it. In the unzipped folder, there should three .dmg files that allow you to download the three Pupil-labs applications.

Using the software

Pupil Labs comes with three different applications, Pupil Capture, Pupil Service, and Pupil Player.



Pupil Service

Pupil service lets you look at each of the eye cameras to make sure they are set up properly for capturing data later.

On launch, you should see a screen that looks similar to the following image:

	Pupil Service		
Settings			
Interface size 1.0	₽		
Reset window size			
Pupil Remote address 134.82.166.26:50020			
Detection & mapping mode 2d	₽		
Detect eye 0	•		
Detect eye 1	•		
Service Version: 1.8.26			
Restart with default settings			

Most of the stuff we actually care about is in the red box. Clicking on the grey circles on the far right will open up the camera feed from the respective eye. From here, we can start calibrating the cameras to get a good image. The window that pops up should look something like the following image:



More information on how to get a good image of your pupil is in the <u>Getting good data</u> section of this guide.

How to navigate this app.

On the far-right side of the screen, there should be four tabs similar to the image below. Hover your mouse over the first tab, which should be labeled "General Setting". Click to open the tab. When opened, the image below should display.

General	\$
Interface Size 0.8 🤤	A
Reset window size	
Flip image display	0
Mode ROI 🤤	C ^{2D}
Click and drag on the blue circles to adjust the region of interest. The region should be as small as possible, but large enough to capture all pupil movements.	
Detection method C++ 2d dete 🤝	

The most useful options here are the Flip Image display button and the Mode dropdown menu. The Flip image display button flips the camera image if it is upside down. The Mode dropdown menu on the other hand, has 3 options: Camera Image, ROI, and Algorithm. Camera Image will show the raw feed from the camera, and Algorithm is helpful when adjusting the filter parameters.

The next tab we will be looking at is the Pupil Detector 2D tab. It should be the bottom tab. Hover over it to confirm you are on the right tab (If your tab says "Pupil Detector 3D", you are still on the right track. If you want to switch to 2d, go to the Pupil Service Menu opened on a different window on your computer and select 2d on the Detection and Mapping Mode dropdown menu). Click on the Pupil Detector 2D tab to open. You should now see the image below:

Pupil Detector 2D	•
Switch to the algorithm display mode to see a visualization of pupil detection parameters overlaid on the eye video. Adjust the pupil intensity range so that the pupil is	
Pupil intensity range 23	(²⁰
Pupil min 68	
Pupil max 150	
Open debug window	

Note: If you do not have a 4th tab option or your tab says Pupil Detector 3D, go back to the main screen (which should still be open in a separate window) and make sure that under Detection & Mapping mode either 2D or 3D is selected.

This screen will either say Pupil Detector 2D or Pupil Detector 3D depending on what detection & mapping mode was selected, but since the menu differences between these two are minimal, we'll just look at 2D here. Here, you can fiddle with the minimum and maximum detectable size of the pupil as well as the contrast. It is extremely advantageous to be in Algorithm mode here, as it allows you to visualize how the pupil detection changes with the settings. Clicking Open Debug menu should display something that looks similar to the screen on the right in the image below:



Note: this will look different if you are in Pupil Detector 3D

The debug screen is helpful for ensuring that the full range of motion of the pupil is able to be detected, and the movements of the pupil in the camera image and the tracked pupil should correlate.

Note: We won't be looking into the other two tabs but feel free to explore them on your own. They are useful for tinkering with the source image.

Once you have used pupil service on both eyes and ensured that both pupils are being tracked cleanly, you can start recording data with Pupil Capture. Refer to Getting Good Data to see if your pupils are tracked cleanly.

Pupil Capture

Upon opening Pupil Capture, you should see something similar to the following image:



It is normal (at least on a Mac computer) for the smaller Pupil Capture windows to pop up, one for each eye. If you accidentally closed one, just click the circle in the General tab to reopen the eye detection.



We start by calibrating. Click the circle with the 'C' on the left hand side of the screen and you should see a marker in the middle of the screen with a red dot. Look at the marker until the dot turns green, and repeat this process until all 5 markers have turned green.



Now that the calibration is completed, the red dot in the middle of the screen should correspond to where in the environment the pupils are focused. If this is not the case, try re-calibrating or clicking the 'T' in the circle to preform an accuracy test.

There are several tabs on the far right of the screen, and 'General' is typically opened on launch. There are 8 icons from top to bottom, and then a horizontal line followed by 2 more icons. Any icon below the line is a plugin, and there is a tab specifically for adding and removing these plugins from the tab bar. I will give a brief overview of the tabs, however most or relatively self-explanatory or can be learned by briefly playing around with their features.

The general tab functions similarly to the tab of the same name in Pupil Service, with the exception of audio-mode allowing the user to change what kind of audio will be played during calibrations and recordings.

General		*
Interface size 1.0		⇒ (#
Reset window size		
audio_mode voice and	d sound	
sound or voice on silent	nly ly	
detection & manning m	node disabled	
detection a mapping m	2d 3d	
Detect eye 0		• •
Detect eye 1		
Capture Version: 1.8.26		•
Restart with default se	ttings	
		S

The next tab is the Plugin Manager, which allows the user to turn on or off icons for specific plugins. There are quite a few options and I encourage playing around with some, but I will go into certain plugins later.



System Graphs dictates what information is displayed on the top left portion of the screen.

System Graphs	
Display CPU usage	
Display frames per second	
Display confidence for eye 0	
Display confidence for eye 1	
Display pupil diameter for eye 0	
Display pupil diameter for eye 1	
Minimum pupil diameter 0.00	
Maximum pupil diameter 8.00	

The next 2 tabs are for image settings, which I rarely used and won't go over since they are pretty intuitive. The Recorder tab allows the user to change the settings for how videos are saved during recordings.

Recorder	
Pupil recording "path_to_recor increasing num name to create	s are saved like this: dings/recording_session_name/nnn" where "nnn" is an aber to avoid overwrites. You can use "/" in your session subdirectories.
Recordings are here but note t	saved to "~/pupil_recordings". You can change the path hat invalid input will be ignored.
Path to record	ings /Users/joaniegannon/recordings
Recording sess	sion name 2018_11_02
Request additi	onal user info
Compression	bigger file, less CPU smaller file, more CPU
Recording the I	raw eye video is optional. We use it for debugging.
Record eye	

The Screen Marker Calibration tab allows for different kinds of calibration that may or may not be easier to preform depending on what setup is being used. Selecting one of these options will give a brief blurb into what they are and how they work. The user can also change calibration settings here.

	ration	•
Calibration Method	Fingertip Calibration HMD Calibration	##
	HMD Calibration 3D	
	Natural Features Calibration	
	Screen Marker Calibration	
	Single Marker Calibration	E
Calibration only cons confidence than the	iders pupil data that has an equal or higher minimum calibration confidence.	0
Minimum calibration	n confidence 0.80	0
Calibrate gaze param	neters using a screen based animation.	0
Monitor Color LCD		0
Use fullscreen		
		۲ س
Marker size 1.00		
Marker size 1.00	•••••••••	

I'm not really sure what purpose this tab has other than debugging, but it looks cool.



As far as the plugins are concerned, most of them are features that may or may not be helpful depending on what the user is attempting to accomplish.

Accuracy Visualizer	
The calibration area (green) is defined as the convex hull of the reference points that were used for calibration. 2D mapping looses accuracy outside of this area. It is recommended to calibrate a big portion of the subject's field of view.	
Visualize mapping error	
The mapping error (orange line) is the angulardistance between mapped pupil positions (red) andtheir corresponding reference points (blue).	
Visualize calibration area	0
Measure gaze mapping accuracy and precision using samples that were collected during calibration. The outlier threshold discards samples with high angular errors.	0
Outlier Threshold [degrees] 5.0	1
Accuracy is calculated as the average angular offset (distance) (in degrees of visual angle) between fixation locations and the corresponding locations of the fixation targets.	
Angular Accuracy 0.347179	°
Precision is calculated as the Root Mean Square (RMS) of the angular distance (in degrees of visual angle) between successive samples during a fixation.	
Angular Precision 0.051177	

When you are ready to record, hit the 'R' button on the left side of the screen. Clicking this button again will stop the recording. Recordings can be viewed using the Pupil Player.



Pupil Player

After using Pupil Capture to record data (by hitting the R button on the left side of the screen), the recorded data should be located at the directory specified (by default: Users/yourUser/recordings)

Opening the app shows the following screen:



Located the saved recording file you wish to play, and drag it onto the window.



The player will show a playback of the recording that looks similar to the following screen:



Click the buttons in the bottom left corner to control the video. The button on the far left will export a .mp4 version of the video, as well as raw data in the form of .csv files. After the export successfully completes, navigate back to the original file, and look in the 'exports' folder within that directory to find the exported data and video.



Getting good data

Good data starts with a good pupil images. To get good pupil readings, make sure the camera is set up such that:

- The full range of motion of the pupil is visible. Make sure that the pupil can still be detected from most if not all angles of vision.
- The eye is in the center of the screen
- There is nothing else in the view of the camera (such as eyebrows).



Good: the full range of motion is visible



Not good: The full range of motion for the pupil is not visible

To get better feedback on weather the eye is properly positioned in the Pupil Service app by using the algorithm view and the debug option. There, parameters such as pupil size and lighting threshold can be adjusted. Additionally, the hardware itself can be adjusted to provide a better view of the pupil. The easiest adjustment is to gently pivot the pupil camera with your fingers. If the camera is not pivoting, you can loosen the screw on the outside of the camera until it is able to move freely while still being tight enough to hold its position while being used.



If these adjustments are still not enough to properly frame the eye in the view of the camera, then the orange extenders can be used.



To attach the extenders, slide the extruding camera out of the side of the glasses, and insert the thin end of the extender into the slot on the camera, like in the image shown below.



Then, slide the extender back onto the glasses where the camera was removed from, as seen in the image below.



Troubleshooting

Here I'll run through some of the issues I've encountered and how I resolved them.

I'm getting frozen images/a grey screen in Pupil Service

Don't worry, you probably didn't break anything, the camera connection probably just came loose. The most likely candidate is the small port that connects the pupil camera to the glasses (seen in the image below), which can sometimes get pulled loose when reorienting the glasses. Gently push the connection back together with your fingertip.



If this did not work, there might be an issue with the USB cable. Otherwise, there might be a loose wire somewhere else on the glasses.

Pupil Capture's calibration is not working

This can be a handful of issues. Firstly, if you are not even able to complete the calibration because the targets don't turn green when looking at them, check to make the pupils are both able to be properly read in Pupil Service. If they are, make sure the head camera is both in focus and is able to see the calibration marks.

If you are able to make it through the calibration, but the calibration still fails, Pupil Capture usually prints out an error message that can be helpful. Make sure both eyes are being detected and are set up properly to be read.

Links and References

Pupil Docs (official documentation): https://docs.pupil-labs.com/#getting-started

Pupil Labs' Github https://github.com/pupil-labs

Pupil Labs' Applications https://github.com/pupil-labs/pupil/releases/latest